

## Nuclear receptor subfamily 2 group C member 2

Catalog No: #AP86169



Package Size: #AP86169-1 50ug #AP86169-2 100ug #AP86169-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear receptor subfamily 2 group C member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Orphan nuclear receptor TR4, Testicular receptor 4, Tr4
Accession No.	Uniprot ID:P55094
Uniprot	P55094
GeneID	50659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Nr2c2

Note: This product is for in vitro research use only