

Mesenteric estrogen-dependent adipogenesis protein

Catalog No: #AP86170



Package Size: #AP86170-1 50ug #AP86170-2 100ug #AP86170-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mesenteric estrogen-dependent adipogenesis protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:A4IFN2
Uniprot	A4IFN2
GeneID	510187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MEDAG

Note: This product is for in vitro research use only