

5-oxoprolinase

Catalog No: #AP86171



Package Size: #AP86171-1 50ug #AP86171-2 100ug #AP86171-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	5-oxoprolinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	5-oxo-L-prolinase,5-OPase,Pyroglutamase
Accession No.	Uniprot ID:Q75WB5
Uniprot	Q75WB5
GeneID	408006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OPLAH

Note: This product is for in vitro research use only