

Matrix Gla protein

Catalog No: #AP86173



Package Size: #AP86173-1 50ug #AP86173-2 100ug #AP86173-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrix Gla protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MGP, Cell growth-inhibiting gene 36 protein, MGLAP, GIG36
Accession No.	Uniprot ID:P08493
Uniprot	P08493
GeneID	4256;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MGP

Note: This product is for in vitro research use only