

Matrix metalloproteinase-14

Catalog No: #AP86175

Package Size: #AP86175-1 50ug #AP86175-2 100ug #AP86175-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MMP-14, Membrane-type matrix metalloproteinase 1, MT-MMP 1, Membrane-type-1 matrix metalloproteinase, MT1-MMP
Accession No.	Uniprot ID:Q9XT90
Uniprot	Q9XT90
GeneID	397471;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MMP14

Note: This product is for in vitro research use only