

Matrix metalloproteinase-21

Catalog No: #AP86176



Package Size: #AP86176-1 50ug #AP86176-2 100ug #AP86176-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MMP-21
Accession No.	Uniprot ID:Q8K3F2
Uniprot	Q8K3F2
GeneID	214766;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mmp21

Note: This product is for in vitro research use only