

Matrix metalloproteinase-20

Catalog No: #AP86177



Package Size: #AP86177-1 50ug #AP86177-2 100ug #AP86177-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Matrix metalloproteinase-20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MMP-20,Enamel metalloproteinase,Enamelysin
Accession No.	Uniprot ID:O18767
Uniprot	O18767
GeneID	281916;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MMP20

Note: This product is for in vitro research use only