

Dual specificity mitogen-activated protein kinase kinase 1

Catalog No: #AP86178

Package Size: #AP86178-1 50ug #AP86178-2 100ug #AP86178-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase kinase 1,ERK activator kinase 1,MAPK/ERK kinase 1,MEK 1,MEK1,PRKMK1
Accession No.	Uniprot ID:Q02750
Uniprot	Q02750
GeneID	5604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MAP2K1

Note: This product is for in vitro research use only