

Pro-neuropeptide Y

Catalog No: #AP86180



Package Size: #AP86180-1 50ug #AP86180-2 100ug #AP86180-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Pro-neuropeptide Y |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Accession No. | Uniprot ID:Q6RUW3 |
| Uniprot | Q6RUW3 |
| GeneID | 504216; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NPY

Note: This product is for in vitro research use only