

Neutrophil gelatinase-associated lipocalin

Catalog No: #AP86181



Package Size: #AP86181-1 50ug #AP86181-2 100ug #AP86181-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neutrophil gelatinase-associated lipocalin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	NGAL,Alpha-2-microglobulin-related protein,Alpha-2U globulin-related protein,Lipocalin-2,Siderocalin LCN2,p25
Accession No.	Uniprot ID:P30152
Uniprot	P30152
GeneID	170496;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lcn2

Note: This product is for in vitro research use only