

Numb-like protein

Catalog No: #AP86182



Package Size: #AP86182-1 50ug #AP86182-2 100ug #AP86182-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Numb-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Numb-related protein, Numb-R
Accession No.	Uniprot ID:Q9Y6R0
Uniprot	Q9Y6R0
GeneID	9253;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NUMBL

Note: This product is for in vitro research use only