

Nitric oxide synthase, endothelial

Catalog No: #AP86184



Package Size: #AP86184-1 50ug #AP86184-2 100ug #AP86184-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nitric oxide synthase, endothelial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Constitutive NOS,cNOS,EC-NOS,Endothelial NOS,eNOS,NOS type III,NOSIII,NOS
Accession No.	Uniprot ID:Q9TUX8
Uniprot	Q9TUX8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NOS3

Note: This product is for in vitro research use only