

## Neurogenic locus notch homolog protein 1

Catalog No: #AP86185



Package Size: #AP86185-1 50ug #AP86185-2 100ug #AP86185-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurogenic locus notch homolog protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Notch 1
Accession No.	Uniprot ID:Q07008
Uniprot	Q07008
GeneID	25496;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Notch1

Note: This product is for in vitro research use only