

Neuronal calcium sensor 1

Catalog No: #AP86190



Package Size: #AP86190-1 50ug #AP86190-2 100ug #AP86190-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuronal calcium sensor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	NCS-1,Frequenin homolog,Frequenin-like protein,Frequenin-like ubiquitous protein,FLUP,FREQ
Accession No.	Uniprot ID:P62167
Uniprot	P62167
GeneID	396336;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NCS1

Note: This product is for in vitro research use only