

NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial

Catalog No: #AP86191

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86191-1 50ug #AP86191-2 100ug #AP86191-3 1mg

Description

Product Name	NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Complex I-23kD, CI-23kD, NADH-ubiquinone oxidoreductase 23 kDa subunit, TYKY subunit
Accession No.	Uniprot ID:O00217
Uniprot	O00217
GeneID	4728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NDUFS8

Note: This product is for in vitro research use only