

Substance-P receptor

Catalog No: #AP86193



Package Size: #AP86193-1 50ug #AP86193-2 100ug #AP86193-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Substance-P receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | SPR,NK-1 receptor,NK-1R,Tachykinin receptor 1,NK1R,TAC1R |
| Accession No. | Uniprot ID:P25103 |
| Uniprot | P25103 |
| GeneID | 6869; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TACR1

Note: This product is for in vitro research use only