

Nuclear receptor subfamily 4 group A member 1

Catalog No: #AP86198



Package Size: #AP86198-1 50ug #AP86198-2 100ug #AP86198-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 4 group A member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Nuclear hormone receptor NUR/77,Nuclear protein N10,Orphan nuclear receptor HMR,Gfrp,Hmr,N10,Nur77
Accession No.	Uniprot ID:P12813
Uniprot	P12813
GeneID	15370;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nr4a1

Note: This product is for in vitro research use only