

Substance-P receptor

Catalog No: #AP86200



Package Size: #AP86200-1 50ug #AP86200-2 100ug #AP86200-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Substance-P receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	SPR,NK-1 receptor,NK-1R,Tachykinin receptor 1
Accession No.	Uniprot ID:Q5DUB3
Uniprot	Q5DUB3
GeneID	403815;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TACR1

Note: This product is for in vitro research use only