

Mu-type opioid receptor

Catalog No: #AP86205



Package Size: #AP86205-1 50ug #AP86205-2 100ug #AP86205-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mu-type opioid receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	M-OR-1,Mu opiate receptor,Mu opioid receptor,MOP,MOR1
Accession No.	Uniprot ID:P35372
Uniprot	P35372
GeneID	4988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OPRM1

Note: This product is for in vitro research use only