

Prefoldin subunit 2

Catalog No: #AP86213

Package Size: #AP86213-1 50ug #AP86213-2 100ug #AP86213-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prefoldin subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pfd2
Accession No.	Uniprot ID:O70591
Uniprot	O70591
GeneID	18637;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pfdn2

Note: This product is for in vitro research use only