

Phosphoprotein associated with glycosphingolipid-enriched microdomains 1

Catalog No: #AP86218

Package Size: #AP86218-1 50ug #AP86218-2 100ug #AP86218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoprotein associated with glycosphingolipid-enriched microdomains 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Csk-binding protein, Transmembrane adapter protein PAG, Transmembrane phosphoprotein Cbp, CBP, PAG
Accession No.	Uniprot ID: Q9NWQ8
Uniprot	Q9NWQ8
GeneID	55824;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PAG1

Note: This product is for in vitro research use only