

## Serine/threonine-protein kinase mTOR

Catalog No: #AP86223



Package Size: #AP86223-1 50ug #AP86223-2 100ug #AP86223-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase mTOR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FK506-binding protein 12-rapamycin complex-associated protein 1,FKBP12-rapamycin complex-associated protein,Mammalian target of rapamycin,mTOR,Mechanistic target of rapamycin,Rapamycin target protein 1,RAPT1,Frap,Frap1
Accession No.	Uniprot ID:Q9JLN9
Uniprot	Q9JLN9
GeneID	56717;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mtor

Note: This product is for in vitro research use only