

Myosin regulatory light polypeptide 9

Catalog No: #AP86224



Package Size: #AP86224-1 50ug #AP86224-2 100ug #AP86224-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin regulatory light polypeptide 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	20 kDa myosin light chain,LC20,MLC-2C,Myosin RLC,Myosin regulatory light chain 2, smooth muscle isoform,Myosin regulatory light chain 9,Myosin regulatory light chain MRLC1,MLC2,MRLC1,MYRL2
Accession No.	Uniprot ID:P24844
Uniprot	P24844
GeneID	10398;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYL9

Note: This product is for in vitro research use only