

P2X purinoceptor 7

Catalog No: #AP86228



Package Size: #AP86228-1 50ug #AP86228-2 100ug #AP86228-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2X purinoceptor 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	P2X7,ATP receptor,P2Z receptor,Purinergic receptor,P2x7
Accession No.	Uniprot ID:Q9Z1M0
Uniprot	Q9Z1M0
GeneID	18439;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P2rx7

Note: This product is for in vitro research use only