

Plasminogen activator inhibitor 2

Catalog No: #AP86230



Package Size: #AP86230-1 50ug #AP86230-2 100ug #AP86230-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Plasminogen activator inhibitor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PAI-2, Monocyte Arg-serpin, Placental plasminogen activator inhibitor, Serpin B2, Urokinase inhibitor, PAI2, PLANH2
Accession No.	Uniprot ID:P05120
Uniprot	P05120
GeneID	5055;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SERPINB2

Note: This product is for in vitro research use only