

## NACHT, LRR and PYD domains-containing protein 3

Catalog No: #AP86233



Package Size: #AP86233-1 50ug #AP86233-2 100ug #AP86233-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NACHT, LRR and PYD domains-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cold autoinflammatory syndrome 1 protein homolog, Cryopyrin, Mast cell maturation-associated-inducible protein 1, PYRIN-containing APAF1-like protein 1, Cias1, Mmig1, Nalp3, Pypaf1
Accession No.	Uniprot ID:Q8R4B8
Uniprot	Q8R4B8
GeneID	216799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Nlrp3

Note: This product is for in vitro research use only