

Vasopressin-neurophysin 2-copeptin

Catalog No: #AP86234



Package Size: #AP86234-1 50ug #AP86234-2 100ug #AP86234-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vasopressin-neurophysin 2-copeptin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AVP-NPII
Accession No.	Uniprot ID:P01180
Uniprot	P01180
GeneID	280728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: AVP

Note: This product is for in vitro research use only