

## Neurabin-2

Catalog No: #AP86236



Package Size: #AP86236-1 50ug #AP86236-2 100ug #AP86236-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurabin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neurabin-II,Protein phosphatase 1 regulatory subunit 9B,Spinophilin,PPP1R6
Accession No.	Uniprot ID:Q96SB3
Uniprot	Q96SB3
GeneID	84687;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPP1R9B

Note: This product is for in vitro research use only