

# Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial

Catalog No: #AP86238

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86238-1 50ug #AP86238-2 100ug #AP86238-3 1mg

## Description

Product Name	Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDHE1-A type I,Pdha-1
Accession No.	Uniprot ID:P35486
Uniprot	P35486
GeneID	18597;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pdha1

Note: This product is for in vitro research use only