

Platelet-derived growth factor D

Catalog No: #AP86239



Package Size: #AP86239-1 50ug #AP86239-2 100ug #AP86239-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PDGF-D,Iris-expressed growth factor,Spinal cord-derived growth factor B,SCDGF-B,legf,Scdgfb
Accession No.	Uniprot ID:Q9EQT1
Uniprot	Q9EQT1
GeneID	66018;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdgfd

Note: This product is for in vitro research use only