

## S-adenosylmethionine synthase isoform type-2

Catalog No: #AP86243



Package Size: #AP86243-1 50ug #AP86243-2 100ug #AP86243-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	S-adenosylmethionine synthase isoform type-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AdoMet synthase 2, Methionine adenosyltransferase 2, MAT 2, Methionine adenosyltransferase II, MAT-II, Ams2
Accession No.	Uniprot ID:P18298
Uniprot	P18298
GeneID	171347;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Mat2a

Note: This product is for in vitro research use only