

Melanoma-derived growth regulatory protein

Catalog No: #AP86245



Package Size: #AP86245-1 50ug #AP86245-2 100ug #AP86245-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Melanoma-derived growth regulatory protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Melanoma inhibitory activity protein
Accession No.	Uniprot ID:Q16674
Uniprot	Q16674
GeneID	8190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MIA

Note: This product is for in vitro research use only