

Orexin receptor type 2

Catalog No: #AP86246



Package Size: #AP86246-1 50ug #AP86246-2 100ug #AP86246-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orexin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Ox-2-R,Hypocretin receptor type 2
Accession No.	Uniprot ID:P56719
Uniprot	P56719
GeneID	25605;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hcrtr2

Note: This product is for in vitro research use only