

Pyrroline-5-carboxylate reductase 2

Catalog No: #AP86247



Package Size: #AP86247-1 50ug #AP86247-2 100ug #AP86247-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyrroline-5-carboxylate reductase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	P5C reductase 2
Accession No.	Uniprot ID:Q6AY23
Uniprot	Q6AY23
GeneID	364064;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pycr2

Note: This product is for in vitro research use only