

Phospholipase A2, membrane associated

Catalog No: #AP86248



Package Size: #AP86248-1 50ug #AP86248-2 100ug #AP86248-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase A2, membrane associated
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GIIC sPLA2, Group IIA phospholipase A2, Non-pancreatic secretory phospholipase A2, NPS-PLA2, Phosphatidylcholine 2-acylhydrolase 2A, PLA2B, PLA2L, RASF-A
Accession No.	Uniprot ID:P14555
Uniprot	P14555
GeneID	5320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PLA2G2A

Note: This product is for in vitro research use only