

Pituitary adenylate cyclase-activating polypeptide type I receptor

Catalog No: #AP86249

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86249-1 50ug #AP86249-2 100ug #AP86249-3 1mg

Description

| | |
|--------------------|--|
| Product Name | Pituitary adenylate cyclase-activating polypeptide type I receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | PACAP type I receptor |
| Accession No. | Uniprot ID:P32215 |
| Uniprot | P32215 |
| GeneID | 24167; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Adcyap1r1

Note: This product is for in vitro research use only