

## Metal regulatory transcription factor 1

Catalog No: #AP86252



Package Size: #AP86252-1 50ug #AP86252-2 100ug #AP86252-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metal regulatory transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MRE-binding transcription factor, Transcription factor MTF-1, Mtf-1
Accession No.	Uniprot ID:Q07243
Uniprot	Q07243
GeneID	17764;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Mtf1

Note: This product is for in vitro research use only