

Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit

Catalog No: #AP86253

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86253-1 50ug #AP86253-2 100ug #AP86253-3 1mg

Description

Product Name	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 48 kDa subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	DDOST 48 kDa subunit,OST48
Accession No.	Uniprot ID:Q29381
Uniprot	Q29381
GeneID	397385;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DDOST

Note: This product is for in vitro research use only