

Neuronal regeneration-related protein

Catalog No: #AP86254



Package Size: #AP86254-1 50ug #AP86254-2 100ug #AP86254-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuronal regeneration-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Neuronal protein 3.1,Protein p311,P311
Accession No.	Uniprot ID:Q80Z34
Uniprot	Q80Z34
GeneID	338475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nrep

Note: This product is for in vitro research use only