

Neurotrophin-4

Catalog No: #AP86255



Package Size: #AP86255-1 50ug #AP86255-2 100ug #AP86255-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NT-4,Neurotrophin-5,NT-5,Neutrophic factor 4,Ntf5
Accession No.	Uniprot ID:Q80VU4
Uniprot	Q80VU4
GeneID	78405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ntf4

Note: This product is for in vitro research use only