

OX-2 membrane glycoprotein

Catalog No: #AP86263



Package Size: #AP86263-1 50ug #AP86263-2 100ug #AP86263-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OX-2 membrane glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD200,MOX1,MOX2,My033
Accession No.	Uniprot ID:P41217
Uniprot	P41217
GeneID	4345;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD200

Note: This product is for in vitro research use only