

## Proteinase-activated receptor 1

Catalog No: #AP86265



Package Size: #AP86265-1 50ug #AP86265-2 100ug #AP86265-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Proteinase-activated receptor 1                        |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Ms   |
| Other Names        | PAR-1,Thrombin receptor,Cf2r,Par1                      |
| Accession No.      | Uniprot ID:P30558                                      |
| Uniprot            | P30558   |
| GeneID             | 14062;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:F2r

Note: This product is for in vitro research use only