

E3 SUMO-protein ligase PIAS3

Catalog No: #AP86267



Package Size: #AP86267-1 50ug #AP86267-2 100ug #AP86267-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 SUMO-protein ligase PIAS3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein inhibitor of activated STAT protein 3
Accession No.	Uniprot ID:Q9Y6X2
Uniprot	Q9Y6X2
GeneID	10401;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PIAS3

Note: This product is for in vitro research use only