

## E3 SUMO-protein ligase PIAS3

Catalog No: #AP86268



Package Size: #AP86268-1 50ug #AP86268-2 100ug #AP86268-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | E3 SUMO-protein ligase PIAS3   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Rt   |
| Other Names        | KChAP,Potassium channel-associated protein,Protein inhibitor of activated STAT protein 3 |
| Accession No.      | Uniprot ID:O70260  |
| Uniprot            | O70260   |
| GeneID             | 83614;   |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pias3

Note: This product is for in vitro research use only