

Numb-like protein

Catalog No: #AP86274



Package Size: #AP86274-1 50ug #AP86274-2 100ug #AP86274-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Numb-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Nbl
Accession No.	Uniprot ID:A1L1I3
Uniprot	A1L1I3
GeneID	292732;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Numb1

Note: This product is for in vitro research use only