

Protein orai-3

Catalog No: #AP86275



Package Size: #AP86275-1 50ug #AP86275-2 100ug #AP86275-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein orai-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transmembrane protein 142C, TMEM142C
Accession No.	Uniprot ID:Q9BRQ5
Uniprot	Q9BRQ5
GeneID	93129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ORAI3

Note: This product is for in vitro research use only