

2'-5'-oligoadenylate synthase 1

Catalog No: #AP86276



Package Size: #AP86276-1 50ug #AP86276-2 100ug #AP86276-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-5'-oligoadenylate synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	(2-5')oligo(A) synthase 1,E18/E16,p46/p42 OAS,OIAS
Accession No.	Uniprot ID:P00973
Uniprot	P00973
GeneID	4938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OAS1

Note: This product is for in vitro research use only