

P2X purinoceptor 7

Catalog No: #AP86281



Package Size: #AP86281-1 50ug #AP86281-2 100ug #AP86281-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2X purinoceptor 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	P2X7,ATP receptor,P2Z receptor,Purinergic receptor
Accession No.	Uniprot ID:Q99572
Uniprot	Q99572
GeneID	5027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P2RX7

Note: This product is for in vitro research use only