

NEDD8 ultimate buster 1

Catalog No: #AP86284



Package Size: #AP86284-1 50ug #AP86284-2 100ug #AP86284-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NEDD8 ultimate buster 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Negative regulator of ubiquitin-like proteins 1,Protein BS4,Nyren18
Accession No.	Uniprot ID:P54729
Uniprot	P54729
GeneID	53312;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Nub1

Note: This product is for in vitro research use only