

Nuclear receptor subfamily 1 group I member 2

Catalog No: #AP86290



Package Size: #AP86290-1 50ug #AP86290-2 100ug #AP86290-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear receptor subfamily 1 group I member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Orphan nuclear receptor PAR1, Orphan nuclear receptor PXR, Pregnane X receptor, Steroid and xenobiotic receptor, SXR, PXR
Accession No.	Uniprot ID: O75469
Uniprot	O75469
GeneID	8856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: NR112

Note: This product is for in vitro research use only